

Human ATR/Serine/threonine-protein kinase ATR ELISA Kit

Catalog No: #EK19714



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ATR/Serine/threonine-protein kinase ATR ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Serine/threonine-protein kinase ATR
Target Name	Serine/threonine-protein kinase ATR
Other Names	Ataxia telangiectasia and Rad3-related protein,FRAP-related protein 1,FRP1
Accession No.	Q13535

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.066ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ATR

Note: This product is for in vitro research use only and is not intended for use in humans or animals.