

Human CD8B/T-cell surface glycoprotein CD8 beta chain ELISA Kit

Catalog No: #EK19722



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CD8B/T-cell surface glycoprotein CD8 beta chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human T-cell surface glycoprotein CD8 beta chain
Target Name	T-cell surface glycoprotein CD8 beta chain
Other Names	CD8b,CD8B1
Accession No.	P10966

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:24pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CD8B

Note: This product is for in vitro research use only and is not intended for use in humans or animals.