

Rat Cnr1/Cannabinoid receptor 1 ELISA Kit

Catalog No: #EK19725



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cnr1/Cannabinoid receptor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cannabinoid receptor 1
Target Name	Cannabinoid receptor 1
Other Names	CB-R,Brain-type cannabinoid receptor,Skr6
Accession No.	P20272

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.089ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cnr1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.