

Mouse Bnip3/BCL2/adenovirus E1B 19 kDa protein-interacting protein 3 ELISA Kit

Catalog No: #EK19729



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Bnip3/BCL2/adenovirus E1B 19 kDa protein-interacting protein 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Target Name	BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Other Names	Nip3
Accession No.	O55003

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.055ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bnip3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.