

# Mouse Akr7a2/Aflatoxin B1 aldehyde reductase member 2 ELISA Kit

Catalog No: #EK19731



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Akr7a2/Aflatoxin B1 aldehyde reductase member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Aflatoxin B1 aldehyde reductase member 2
Target Name	Aflatoxin B1 aldehyde reductase member 2
Other Names	Succinic semialdehyde reductase, SSA reductase, Afar, Akr7a5
Accession No.	Q8CG76

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Akr7a2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.