

# Human AQP12A/Aquaporin-12A ELISA Kit

Catalog No: #EK19739



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human AQP12A/Aquaporin-12A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aquaporin-12A
Target Name	Aquaporin-12A
Other Names	AQP-12,AQP12,AQPX2
Accession No.	Q8IXF9

## Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/ml Sensitivity:0.36ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:AQP12A

Note: This product is for in vitro research use only and is not intended for use in humans or animals.