

# Rat Cth/Cystathionine gamma-lyase ELISA Kit

Catalog No: #EK19741



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cth/Cystathionine gamma-lyase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cystathionine gamma-lyase
Target Name	Cystathionine gamma-lyase
Other Names	Cysteine-protein sulfhydrase,Gamma-cystathionase,Probasin-related antigen,PRB-RA
Accession No.	P18757

## Application Details

Assay Type: Sandwich Detect Range:3.12-200ng/mL Sensitivity:1.62ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cth

Note: This product is for in vitro research use only and is not intended for use in humans or animals.