

# Mouse Apcs/Serum amyloid P-component ELISA Kit

Catalog No: #EK19744



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Apcs/Serum amyloid P-component ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Serum amyloid P-component
Target Name	Serum amyloid P-component
Other Names	SAP,Ptx2,Sap
Accession No.	P12246

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 40.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Apcs

Note: This product is for in vitro research use only and is not intended for use in humans or animals.