

Dog CCL17/C-C motif chemokine 17 ELISA Kit

Catalog No: #EK19748



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Dog CCL17/C-C motif chemokine 17 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Dog |
| Specificity | Natural and recombinant Dog C-C motif chemokine 17 |
| Target Name | C-C motif chemokine 17 |
| Other Names | CC chemokine TARC, Small-inducible cytokine A17, Thymus and activation-regulated chemokine, TARC |
| Accession No. | Q95N01 |

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 40.7pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL17

Note: This product is for in vitro research use only and is not intended for use in humans or animals.