

BCL-2 (Phospho-Thr74) Antibody

Catalog No: #12745

Package Size: #12745-1 50ul #12745-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BCL-2 (Phospho-Thr74) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	BCL-2 (Phospho-Thr74) Antibody detects endogenous levels of BCL-2 only when phosphorylated at Thr74
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human BCL-2 (Phospho-Thr74)
Target Name	BCL-2
Modification	Phospho
Other Names	BCL2, Apoptosis regulator Bcl-2, B-cell CLL/lymphoma 2, B-cell lymphoma protein 2, Bcl-2, PPP1R50
Accession No.	Swiss-Prot#: P10415NCBI Gene ID: 596
Target Species	human
Calculated MW	28kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.