

# Human CCL23/C-C motif chemokine 23 ELISA Kit

Catalog No: #EK19750



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CCL23/C-C motif chemokine 23 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C motif chemokine 23
Target Name	C-C motif chemokine 23
Other Names	CK-beta-8,CKB-8,Macrophage inflammatory protein 3,MIP-3,Myeloid progenitor inhibitory factor 1,MPIF-1,Small-inducible cytokine A23,MIP3,MPIF1,SCYA23
Accession No.	P55773

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/ML Sensitivity: 15pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CCL23

Note: This product is for in vitro research use only and is not intended for use in humans or animals.