

Human CFC1/Cryptic protein ELISA Kit

Catalog No: #EK19757



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CFC1/Cryptic protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cryptic protein
Target Name	Cryptic protein
Other Names	Cryptic family protein 1
Accession No.	P0CG37

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 41.1pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CFC1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.