

# Dog CASP1/Caspase-1 ELISA Kit

Catalog No: #EK19761



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog CASP1/Caspase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Caspase-1
Target Name	Caspase-1
Other Names	CASP-1, Interleukin-1 beta convertase, IL-1BC, Interleukin-1 beta-converting enzyme, ICE, IL-1 beta-converting enzyme, p45, IL1BC
Accession No.	Q9MZV7

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.75pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CASP1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.