

# Mouse Clec4e/C-type lectin domain family 4 member E ELISA Kit

Catalog No: #EK19767



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mouse Clec4e/C-type lectin domain family 4 member E ELISA Kit                        |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Ms   |
| Specificity        | Natural and recombinant Mouse C-type lectin domain family 4 member E                 |
| Target Name        | C-type lectin domain family 4 member E   |
| Other Names        | C-type lectin superfamily member 9,Macrophage-inducible C-type lectin,Clecsf9,Mincle |
| Accession No.      | Q9R0Q8   |

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.08ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Clec4e

Note: This product is for in vitro research use only and is not intended for use in humans or animals.