

Rat Clec4e/C-type lectin domain family 4 member E ELISA Kit

Catalog No: #EK19768



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Clec4e/C-type lectin domain family 4 member E ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat C-type lectin domain family 4 member E
Target Name	C-type lectin domain family 4 member E
Other Names	C-type lectin superfamily member 9, Macrophage-inducible C-type lectin, Clecsf9, Mincle
Accession No.	Q67EQ1

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.081ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Clec4e

Note: This product is for in vitro research use only and is not intended for use in humans or animals.