

Human ANK2/Ankyrin-2 ELISA Kit

Catalog No: #EK19773



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ANK2/Ankyrin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Ankyrin-2
Target Name	Ankyrin-2
Other Names	ANK-2,Ankyrin-B,Brain ankyrin,Non-erythroid ankyrin
Accession No.	Q01484

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.41ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ANK2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.