

Human CLCF1/Cardiotrophin-like cytokine factor 1 ELISA Kit

Catalog No: #EK19778



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Human CLCF1/Cardiotrophin-like cytokine factor 1 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Cardiotrophin-like cytokine factor 1 |
| Target Name | Cardiotrophin-like cytokine factor 1 |
| Other Names | B-cell-stimulating factor 3,BSF-3,Novel neurotrophin-1,NNT-1,BSF3,CLC,NNT1 |
| Accession No. | Q9UBD9 |

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CLCF1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.