

PNKP (Phospho-Ser114) Antibody

Catalog No: #12775

Package Size: #12775-1 50ul #12775-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | PNKP (Phospho-Ser114) Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide. |
| Applications | WB |
| Species Reactivity | Hu |
| Specificity | PNKP (Phospho-Ser114) Antibody detects endogenous levels of PNKP only when phosphorylated at Ser114 |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human PNKP (Phospho-Ser114) |
| Target Name | PNKP |
| Modification | Phospho |
| Other Names | PNKP, EIEE10, DNA 5-kinase/3-phosphatase, PNK, MCSZ |
| Accession No. | Swiss-Prot#: Q96T60NCBI Gene ID: 11284 |
| Target Species | human |
| Calculated MW | 60kd |
| Concentration | 1.0mg/ml |
| Formulation | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at -20°C |

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.