

SYN2 (Phospho-Ser425) Antibody

Catalog No: #12784



Package Size: #12784-1 50ul #12784-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SYN2 (Phospho-Ser425) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SYN2 (Phospho-Ser425) Antibody detects endogenous levels of SYN2 only when phosphorylated at Ser425
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SYN2 (Phospho-Ser425)
Target Name	SYN2
Modification	Phospho
Other Names	SYN2, Synapsin II, SYNII, SYNIIb, Synapsin-2, SYNIIa
Accession No.	Swiss-Prot#: Q92777NCBI Gene ID: 6854
Target Species	human
Calculated MW	78kd?
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.