

PGK1 (Acetyl-Lys388) Antibody

Catalog No: #11599



Package Size: #11599-1 50ul #11599-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PGK1 (Acetyl-Lys388) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB IHC IP
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous level of PGK1 only when acetylated at Lysine 388.
Immunogen Type	Peptide
Immunogen Description	Peptide sequence around acetylation site of lysine 388(E-D-K(Acetyl)-V-S) derived from Human PGK1.
Other Names	MGC117307,MGC142128,MGC8947,MIG10,PGKA
Accession No.	Swiss-Prot#: P00558 NCBI Gene ID: 5230 NCBI mRNA#: NM_000291.3 NCBI Protein#: NP_000282.1
Calculated MW	45KD
Concentration	1.0mg/mL
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB dilution: 1:500~1:1000

IHC dilution: 1:50~1:100

IP dilution: 1:20~1:50

Product Description

Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

Background

The PGK1 gene encodes phosphoglycerate kinase-1, also known as ATP:3-phosphoglycerate 1-phosphotransferase (EC 2.7.2.3), which catalyzes the reversible conversion of 1,3-diphosphoglycerate to 3-phosphoglycerate during glycolysis, generating one molecule of ATP. It belongs to the phosphoglycerate kinase family and defects in PGK1 are the cause of phosphoglycerate kinase 1 deficiency (PGK1D).

References

Xu Qian, et al. Phosphoglycerate kinase 1 phosphorylates Beclin1 to induce autophagy. Mol Cell. 2017 March 02; 65(5): 917–931.e6. doi:10.1016/j.molcel.2017.01.027. Autho

Published Papers

et al., Identification of a novel non-ATP-competitive protein kinase inhibitor of PGK1 from marine nature products. In *Biochem Pharmacol* on 2021 Jan by Yuying Wang, Lulu Sun,

et al..PMID:33212041, , (2021)

[PMID:33212041](#)

et al., Phosphoglycerate Kinase 1 Phosphorylates Beclin1 to Induce Autophagy.In *Mol Cell* on 2017 Mar 2 by Xu Qian , Xinjian Li,et al..PMID:

28238651, , (2017)

[PMID:28238651](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.