

## MERTK (phospho-Tyr753/685) Antibody

Catalog No: #13320

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MERTK (phospho-Tyr753/685) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of MERTK (phospho-Tyr753/Y685) protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human MERTK.
Target Name	MERTK; TYRO3
Other Names	MERTK; MER; Tyrosine-protein kinase Mer; Proto-oncogene c-Mer; Receptor tyrosine kinase MerTK; TYRO3; BYK; DTK; RSE; SKY; Tyrosine-protein kinase receptor TYRO3; Tyrosine-protein kinase DTK; Tyrosine-protein kinase RSE; Tyrosine-protein kinase SKY; Tyrosine-protein kinase byk
Accession No.	Swiss-Prot#:Q12866; Q06418NCBI Gene ID:10461; 7301
Calculated MW	110KD
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting:1:500 - 1:1000

## Images



Western blot analysis of MERTK (phospho-Tyr753/685) expression in HEK293T, H1688 whole cell lysates.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.