

MRGX3 Antibody

Catalog No: #47834



Package Size: #47834-1 50ul #47834-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MRGX3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of MRGX3 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human MRGX3.
Target Name	MRGPRX3
Other Names	MRGX3; SNSR1; SNSR2; Mas-related G-protein coupled receptor member X3; Sensory neuron-specific G-protein coupled receptor 1/2
Accession No.	Swiss-Prot#:Q96LB0NCBI Gene ID:117195
Calculated MW	36KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of MRGX3 expression in human brain
(A) whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.