

## PABPC5 Antibody

Catalog No: #47844



Package Size: #47844-1 50ul #47844-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PABPC5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB, IHC
Species Reactivity	Hu,Rt
Specificity	Recognizes endogenous levels of PABPC5 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PABPC5.
Target Name	PABPC5
Other Names	PABP5; Polyadenylate-binding protein 5; PABP-5; Poly(A)-binding protein 5
Accession No.	Swiss-Prot#:Q96DU9NCBI Gene ID:140886
Calculated MW	43,39KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

## Application Details


Western blotting:1:500 - 1:1000

Immunohistochemistry:1:50 - 1:100

## Images



Western blot analysis of PABPC5 expression in THP1 (A), PC3 (B), HeLa (C), rat brain (D) whole cell lysates.



Immunohistochemical analysis of PABPC5 staining in human lung cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.