

ZNT4 Antibody

Catalog No: #47855



Package Size: #47855-1 50ul #47855-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNT4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of ZNT4 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human SLC30A4.
Target Name	SLC30A4
Other Names	ZNT4; Zinc transporter 4; ZnT-4; Solute carrier family 30 member 4
Accession No.	Swiss-Prot#:O14863NCBI Gene ID:7782
Calculated MW	47KD
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting:1:500 - 1:1000Immunohistochemistry:1:50 - 1:100

Images



Western blot analysis of RPS27L expression in H1792 whole cell lysates.



Immunohistochemical analysis of ZNT4 staining in human lymph node formalin fixed paraffin embedded tissue section. The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.