

GPR172A Antibody

Catalog No: #47856



Package Size: #47856-1 50ul #47856-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPR172A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of GPR172A protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human GPR172A.
Target Name	SLC52A2
Other Names	GPR172A; PAR1; RFT3; Solute carrier family 52, riboflavin transporter, member 2; Porcine endogenous retrovirus A receptor 1; PERV-A receptor 1; Protein GPR172A; Riboflavin transporter 3; hRFT3
Accession No.	Swiss-Prot#:Q9HAB3NCBI Gene ID:79581
Calculated MW	46KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting: 1:500 - 1:1000

Images



Western blot analysis of GPR172A expression in HEK293T (A), MCF7 (B) whole cell lysates.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.