

TMEM145 Antibody

Catalog No: #47861



Package Size: #47861-1 50ul #47861-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMEM145 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB, IHC
Species Reactivity	Hu,Ms
Specificity	Recognizes endogenous levels of TMEM145 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human TMEM145.
Target Name	TMEM145
Other Names	Transmembrane protein 145
Accession No.	Swiss-Prot#:Q8NBT3NCBI Gene ID:284339
Calculated MW	56KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.


Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of TMEM145 expression in human pancreas (A) whole cell lysates



Immunohistochemical analysis of TMEM145 staining in human brain formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.