

## CD32c Antibody

Catalog No: #47869

Package Size: #47869-1 50ul #47869-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CD32c Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of CD32c protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human CD32c.
Target Name	FCGR2C
Other Names	CD32; FCG2; IGFR2; Low affinity immunoglobulin gamma Fc region receptor II-c; IgG Fc receptor II-c; CDw32; Fc-gamma RII-c; Fc-gamma-RIIc; FcRII-c; CD32
Accession No.	Swiss-Prot#:P31995NCBI Gene ID:9103
Calculated MW	36KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

## Application Details

Western blotting:1:500 - 1:1000

## Images



Western blot analysis of CD32c expression in HEK293T (A), mouse lung (B) whole cell lysates.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.