

CNGB1 Antibody

Catalog No: #47874

Package Size: #47874-1 50ul #47874-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CNGB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of CNGB1 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human CNGB1.
Target Name	CNGB1
Other Names	CNCG2; CNCG3L; CNCG4; RCNC2; Cyclic nucleotide-gated cation channel beta-1; Cyclic nucleotide-gated cation channel 4; CNG channel 4; CNG-4; CNG4; Cyclic nucleotide-gated cation channel gamma; Cyclic nucleotide-gated cation channel modulatory subunit; Cyclic nucleotide-gated channel beta-1; CNG channel beta-1; Glutamic acid-rich protein; GARP
Accession No.	Swiss-Prot#:Q14028NCBI Gene ID:1258
Calculated MW	102KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



"Western blot analysis of CNGB1 expression in CT26 (A), HCT116 (B), A375 (C), LO2 (D), MCF7 (E), K562 (F) whole cell lysates."

Note: This product is for in vitro research use only and is not intended for use in humans or animals.