

GALR2 Antibody

Catalog No: #47876

Package Size: #47876-1 50ul #47876-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GALR2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of GALR2 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human GALR2.
Target Name	GALR2
Other Names	GALNR2; Galanin receptor type 2; GAL2-R; GALR-2
Accession No.	Swiss-Prot#:O43603NCBI Gene ID:8811
Calculated MW	42KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



"Western blot analysis of GALR2 expression in Hela (A), HepG2 (B), mouse brain (C), mouse kidney (D), rat kidney (E) whole cell lysates"

Note: This product is for in vitro research use only and is not intended for use in humans or animals.