

## ADAMTS1 Antibody

Catalog No: #47877



Package Size: #47877-1 50ul #47877-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ADAMTS1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of ADAMTS1 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ADAMTS1.
Target Name	ADAMTS1
Other Names	KIAA1346; METH1; A disintegrin and metalloproteinase with thrombospondin motifs 1; ADAM-TS 1; ADAM-TS1; ADAMTS-1; METH-1
Accession No.	Swiss-Prot#:Q9UHI8NCBI Gene ID:9510
Calculated MW	110KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

## Application Details

Western blotting: 1:500 - 1:1000

## Images



"Western blot analysis of ADAMTS1 expression in A549 (A), H446 (B) whole cell lysates."

Note: This product is for in vitro research use only and is not intended for use in humans or animals.