

MRPS33 Antibody

Catalog No: #47898



Package Size: #47898-1 50ul #47898-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	MRPS33 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms
Specificity	Recognizes endogenous levels of MRPS33 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human MRPS33.
Target Name	MRPS33
Other Names	28S ribosomal protein S33, mitochondrial; MRP-S33; S33mt
Accession No.	Swiss-Prot#:Q9Y291NCBI Gene ID:51650
Calculated MW	12KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Immunohistochemistry:1:50 - 1:100

Images



Western blot analysis of MRPS33 expression in HEK293T (A), Raw264.7 (B) whole cell lysates.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.