

# Kappa Opioid Receptor Antibody

Catalog No: #47902

Package Size: #47902-1 50ul #47902-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Kappa Opioid Receptor Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of Kappa Opioid Receptor protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Kappa Opioid Receptor.
Target Name	OPRK1
Other Names	OPRK1; OPRK; Kappa-type opioid receptor; K-OR-1; KOR-1
Accession No.	Swiss-Prot#:P41145NCBI Gene ID:4986
Calculated MW	42KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

## Application Details

Western blotting:1:500 - 1:1000

## Images



Western blot analysis of Kappa Opioid Receptor expression in HeLa (A), HepG2 (B) whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.