

## ABHD7 Antibody

Catalog No: #47917



Package Size: #47917-1 50ul #47917-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ABHD7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of ABHD7 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ABHD7.
Target Name	EPHX4
Other Names	ABHD7; EPHXRP; Epoxide hydrolase 4; Abhydrolase domain-containing protein 7; Epoxide hydrolase-related protein
Accession No.	Swiss-Prot#:Q8IUS5NCBI Gene ID:253152
Calculated MW	42KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

## Application Details

Western blotting: 1:500 - 1:1000

## Images



Western blot analysis of ABHD7 expression in HEK293T (A), HepG2 (B) whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.