

GABRA6 Antibody

Catalog No: #47918

Package Size: #47918-1 50ul #47918-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GABRA6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of GABRA6 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human GABRA6.
Target Name	GABRA6
Other Names	Gamma-aminobutyric acid receptor subunit alpha-6; GABA(A) receptor subunit alpha-6
Accession No.	Swiss-Prot#:Q16445NCBI Gene ID:2559
Calculated MW	51KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of GABRA6 expression in Jurkat (A), mouse brain (B) whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.