

Rabenosyn 5 Antibody

Catalog No: #47935



Package Size: #47935-1 50ul #47935-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Rabenosyn 5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of Rabenosyn 5 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Rabenosyn 5.
Target Name	ZFYVE20
Other Names	Rabenosyn-5; 110 kDa protein; FYVE finger-containing Rab5 effector protein rabenosyn-5; Zinc finger FYVE domain-containing protein 20
Accession No.	Swiss-Prot#:Q9H1K0NCBI Gene ID:64145
Calculated MW	89KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of Rabenosyn 5 expression in A549 (A), Hela (B), BV2 (C), PMVEC (D) whole cell lysates.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.