

VPS4B Antibody

Catalog No: #47946



Package Size: #47946-1 50ul #47946-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VPS4B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of VPS4B protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human VPS4B.
Target Name	VPS4B
Other Names	SKD1; VPS42; Vacuolar protein sorting-associated protein 4B; Cell migration-inducing gene 1 protein; Suppressor of K(+) transport growth defect 1; Protein SKD1
Accession No.	Swiss-Prot#:O75351NCBI Gene ID:9525
Calculated MW	49KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of VPS4B expression in HEK293T (A), A549 (B), H1688 (C) whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.