

FNBP3 Antibody

Catalog No: #47954

Package Size: #47954-1 50ul #47954-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | FNBP3 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was purified by immunogen affinity chromatography. |
| Applications | WB |
| Species Reactivity | Hu,Ms,Rt |
| Specificity | Recognizes endogenous levels of FNBP3 protein. |
| Immunogen Description | KLH-conjugated synthetic peptide encompassing a sequence within the center region of human FNBP3. |
| Target Name | PRPF40A |
| Other Names | FBP11; FLAF1; FNBP3; HIP10; HYPA; Pre-mRNA-processing factor 40 homolog A; Fas ligand-associated factor 1; Formin-binding protein 11; Formin-binding protein 3; Huntingtin yeast partner A; Huntingtin-interacting protein 10; HIP-10; Huntingtin-interacting protein A; Renal carcinoma antigen NY-REN-6 |
| Accession No. | Swiss-Prot#:O75400NCBI Gene ID:55660 |
| Calculated MW | 108KD |
| Concentration | 1 mg/ml |
| Formulation | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of FNBP3 expression in HeLa (A), mouse brain (B), rat brain (C) whole cell lysates.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.