

HAPLN4 Antibody

Catalog No: #47956

Package Size: #47956-1 50ul #47956-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HAPLN4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of HAPLN4 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human HAPLN4.
Target Name	HAPLN4
Other Names	BRAL2; KIAA1926; Hyaluronan and proteoglycan link protein 4; Brain link protein 2
Accession No.	Swiss-Prot#:Q86UW8NCBI Gene ID:404037
Calculated MW	42KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of HAPLN4 expression in mouse muscle (A), rat kidney (B) whole cell lysates.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.