

TBX2/3 Antibody

Catalog No: #47959

Package Size: #47959-1 50ul #47959-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TBX2/3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of TBX2/3 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human TBX2/3.
Target Name	TBX2; TBX3
Other Names	TBX3; T-box transcription factor TBX3; T-box protein 3; TBX2; T-box transcription factor TBX2; T-box protein 2
Accession No.	Swiss-Prot#:O15119; Q13207NCBI Gene ID:6926
Calculated MW	79KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of TBX2/3 expression in HEK293T (A), Hela (B), H446 (C), rat kidney (D) whole cell lysates.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.