

RNF111 Antibody

Catalog No: #47969

Package Size: #47969-1 50ul #47969-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RNF111 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of RNF111 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human RNF111.
Target Name	RNF111
Other Names	E3 ubiquitin-protein ligase Arkadia; RING finger protein 111
Accession No.	Swiss-Prot#:Q6ZNA4NCBI Gene ID:54778
Calculated MW	109KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

Western blotting:1:500 - 1:1000

Images



Western blot analysis of RNF111 expression in Hela (A), A549 (B), MCF7 (C), CT26 (D), PC12 (E) whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.