

# Goat Anti-Human IgG Secondary Antibody AP Conjugated

Catalog No: #L3023



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #L3023-1 100ul #L3023-2 1.0 ml

## Description

Product Name	Goat Anti-Human IgG Secondary Antibody AP Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	IgG (H+L)
Conjugates	Alkaline Phosphatase
Formulation	0.01M Tris-HCl, 0.25M NaCl, pH 8.0, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

## Application Details

Antibody Format: Whole IgG

ELISA: 1:5,000-1:50,000

Western Blot: 1:5,000-1:50,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.