

Goat Anti-Human IgG Secondary Antibody PE Conjugated

Catalog No: #L30217



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #L30217-1 100ul #L30217-2 1.0 ml

Description

Product Name	Goat Anti-Human IgG Secondary Antibody PE Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	IgG (H+L)
Conjugates	R-Phycoerythrin
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide,
Storage	Store at -20 °C/1 year

Application Details

Antibody Format:F(ab') Fragment

most applications: 1:100 - 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.