

Goat Anti-Mouse IgG Secondary Antibody AP Conjugated

Catalog No: #L3033



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #L3033-1 100ul #L3033-2 1.0 ml

Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody AP Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	IgG (H+L)
Conjugates	Alkaline Phosphatase
Formulation	0.01M Tris-HCl, 0.25M NaCl, pH 8.0, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

ELISA: 1:5,000-1:50,000

Western Blot: 1:5,000-1:50,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.