

Goat Anti-Rabbit IgG Secondary Antibody AF594 Conjugated

Catalog No: #L3017

Package Size: #L3017-1 100ul #L3017-2 1.0 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody AF594 Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	AF594
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:800

Published Papers

el at., TRPM7-Mediated Ca²⁺ Regulates Mussel Settlement through the CaMKK β -AMPK-SGF1 Pathway In Int J Mol Sci On 2023 Mar 11 by Jian He?1, Peng Wang et al.. PMID:36982474, (2023)

[PMID:36982474](#)

el at., MicroRNA-22 Inhibits the Apoptosis of Vascular Smooth Muscle Cell by Targeting p38MAPK α in Vascular Remodeling of Aortic Dissection. In Mol Ther Nucleic Acids on 2020 Aug 25 by Yu Xiao, Yudong Sun, et al.. PMID:33294292, (2020)

[PMID:33294292](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.