

Goat Anti-Rabbit IgG Secondary Antibody APC Conjugated

Catalog No: #L30116



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #L30116-1 100ul #L30116-2 0.5 ml

Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody APC Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	Allophycocyanin (APC)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide,
Storage	Store at -20 °C/1 year

Application Details

Antibody Format:F(ab') Fragment

most applications: 1:100 - 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.