

# Goat Anti-Rabbit IgG Secondary Antibody FITC Conjugated

Catalog No: #L30113

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #L30113-1 100ul #L30113-2 1.0 ml

## Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody FITC Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	Fluorescein (FITC)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

## Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:200

## Published Papers

el at., Clock-modified mesenchymal stromal cells therapy rescues molecular circadian oscillation and age-related bone loss via miR142-3p/Bmal1/YAP signaling axis. In *Cell Death Discov.* 2022 Mar 12 by Sa Cha 1, Jiangyue Wang, et al.. PMID:35279674, , (2022)

[PMID:35279674](#)

el at., Hemophagocytic lymphohistiocytosis as an onset of diffuse large B-cell lymphoma: A case report. In *Oncol Lett* on 2022 Jul 5 by Yueqing Cao, Lang Zou, et al.. PMID: 35949601, , (2022)

[PMID:35949601](#)

el at., Oral *Treponema denticola* Infection Induces A $\epsilon$ <sup>o</sup>Y1 $\uparrow$  ?0 and A $\epsilon$ <sup>o</sup>Y1 $\uparrow$  ?2 Accumulation in the Hippocampus of C57BL/6 Mice. In *J Mol Neurosci* on 2021 Jul by Xinyi Su 1, Zhiqun Tang, et al.. PMID: 33763842, , (2021)

[PMID:33763842](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.