

# Goat Anti-Rabbit IgG Secondary Antibody Rhodamine (RRX) Conjugated

Catalog No: #L30115



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #L30115-1 100ul #L30115-2 1.0 ml

## Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody Rhodamine (RRX) Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	Rhodamine Red-X (RRX)
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

## Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.