

ZNF262 Antibody

Catalog No: #44437



Package Size: #44437-1 50ul #44437-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF262 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms
Specificity	ZNF262 Antibody detects endogenous levels of total ZNF262
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZNF262
Target Name	ZNF262
Other Names	ZMYM4, CDIR, Cell death inhibiting RNA, MYM, Zinc finger MYM-type protein 4, Zinc finger protein 262, ZNF198L3, Zinc finger, MYM-type 4, ZNF262, KIAA0425
Accession No.	Swiss-Prot#: Q5VZL5NCBI Gene ID:9202
Calculated MW	175kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western blot analysis ZNF262 using HepG2 whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.