

ZNF785 Antibody

Catalog No: #44438



Package Size: #44438-1 50ul #44438-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF785 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ZNF785 Antibody detects endogenous levels of total ZNF785
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZNF785
Target Name	ZNF785
Other Names	ZNF785, Zinc finger protein 785
Accession No.	Swiss-Prot#: A8K8V0NCBI Gene ID:146540
Calculated MW	87kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western blot analysis ZNF785 using HepG2 whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.