

HSP70 antibody

Catalog No: #23022

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HSP70 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 308 and 579 of HSP70
Target Name	HSP70
Accession No.	NCBI Gene ID: 3303NCBI mRNA#: NM_005345NCBI Protein#: NP_005336
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 70kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images



Sample (30 ug of whole cell lysate)


A: H1299

B: HeLa

C: Hep G2

7.5% SDS PAGE

Primary antibody diluted at 1: 20000



Confocal immunofluorescence analysis (Olympus FV10i) of paraformaldehyde-fixed A431, using HSP70 1A antibody (Green) at 1: 500 dilution and alpha-tubulin antibody (Red) at 1: 2000.

Background

This intronless gene encodes a 70kDa heat shock protein which is a member of the heat shock protein 70 family. In conjunction with other heat shock proteins, this protein stabilizes existing proteins against aggregation and mediates the folding of newly translated proteins in the cytosol and in organelles. It is also involved in the ubiquitin-proteasome pathway through interaction with the AU-rich element RNA-binding protein 1. The gene is located in the major histocompatibility complex class III region, in a cluster with two closely related genes which encode similar proteins. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.